



## NTM PAGER WITH VORN AND DISPLAY

BOGGLE	
	Article number: BOGGLE
	Size: 62 x 31 x 11 mm
	Function: Pager with audio and visual notifying. VORN intercom. Personal guard.
	Standards: EN300-220-1 EN300-220-2 EN300-220-3 IEEE 802.15.4 EN54-25 IIC
<b>Specifications</b>	
<b>Functional diagram</b>	
<b>Description</b>	<p>The Boggle is enhanced NTM transceiver with VORN capability (audio intercom) and equipped with a shock and temperature sensor in the form of a wristwatch.</p> <p>It can be used as a pager in a nurse call system based on the Ninthway high secure radio network</p> <p>It shows calls from NCA or NCB call stations, and with a simple push sets up a VORN (voice over radio network) communication link with the activated NCA station.</p> <p>The Boggle can also be used as a personal watch station or emergency button with built in intercom facility and movement sensor.</p> <p>It has a real time clock, with multiple alarms that generates information on the display and audio messages.</p> <p>A user can activate an alarm broadcast by ticking on the device.</p> <p>The device carries a shock sensor that will detect free fall or abrupt motions</p>

## NTM PAGER WITH VORN AND DISPLAY

<b>BOGGLE</b>	
	<p>indicating the user is falling.</p> <p>The multi alarm feature provides a set of unique and repeating alarms to remind the user on take some action (taking medicine).</p> <p>The Boggle carries a temperature sensor that can give an estimate of the body temperature. Alarms can be set using the sensor.</p> <p>Alarms, display signals and audio signals can be set using remote programming from a base station.</p>
<b>VORN</b>	<p>Electret microphone with built-in AGC.</p> <p>Ceramic speaker.</p>
<b>Sensors</b>	<p>Temperature sensor MCP9800. -10 to +85 °C.</p> <p>3D Accelerometer ADXL345.</p>
<b>Programming</b>	<p>See application note 8; Manual for the Boggle</p>
<b>Indicators</b>	<p>LD1 Green LED</p> <p>LD2 Yellow LED</p> <p>LD3 Red LED</p> <p>2 x 40 character matrix display</p>
<b>Power supply</b>	<p>CR2450 LiSoCl2 in battery holder at the clip.</p>
<b>Radio parameters</b>	<p>See datasheet en manual NCA</p>
<b>NTM mode</b>	<p>The NCA-mode is activated using function 7.</p>
<b>Mounting instructions</b>	<p>Preferably wear on left wrist.</p>
<b>Additional information</b>	<p>Datasheet NTM_3</p> <p>Application note 1; Programming the NTM.</p> <p>Application note 2; Ninthway high secure radio network.</p> <p>Application note 8: Manual for the boggle.</p>