


NURSE CALL BINARY STATION


NCB



Article number:	SPI+NCB_1
Size:	60 x 42 x 18 mm
Function:	Interface between on/off inputs and/or outputs and the Ninthway radio network
Standards:	EN300-220-1 EN300-220-2 EN300-220-3 IEEE 802.15.4 EN54-25 IIC

Specifications

Functional diagram



NURSE CALL BINARY STATION

NCB																	
Description	<p>The NCB_1 is designed as a multipurpose warn station. It uses an 8 channel I2C extender of which four channels are dedicated as switch inputs and four are available as inputs or outputs to be configured at will. Two input switches are already built in.</p> <p>It has two open collector outputs for indicators. Either the built in LED's can be used or external ones.</p> <p>Change of one or more inputs, triggers the NCB, depending on the set flavour, to set the indicators and transmit a frame with a two bytes payload.</p> <ul style="list-style-type: none"> - Byte 0 contains the code for the flavour - Byte 1 contains the status of the extender inputs <p>The NCB comes with a number of pre settable operation modes called flavours that couple the input signals to the indicators.</p> <p>Flavour:</p> <table style="width: 100%; border: none;"> <tr> <td style="width: 50%;">1. Offuse</td> <td>simple I/O device</td> </tr> <tr> <td>2. Door</td> <td>door/window contact</td> </tr> <tr> <td>3. Indicator</td> <td>wireless side indicator</td> </tr> <tr> <td>4. Alarm</td> <td>two button warning station</td> </tr> <tr> <td>5. Nsalarm</td> <td>two button warning station</td> </tr> <tr> <td>6. Triple pull</td> <td>pull cord warning station</td> </tr> <tr> <td>7. External</td> <td>use of external IIC device</td> </tr> <tr> <td>8. Syncalarm</td> <td>synchronized two button warning station to be used with wireless Indicator</td> </tr> </table>	1. Offuse	simple I/O device	2. Door	door/window contact	3. Indicator	wireless side indicator	4. Alarm	two button warning station	5. Nsalarm	two button warning station	6. Triple pull	pull cord warning station	7. External	use of external IIC device	8. Syncalarm	synchronized two button warning station to be used with wireless Indicator
1. Offuse	simple I/O device																
2. Door	door/window contact																
3. Indicator	wireless side indicator																
4. Alarm	two button warning station																
5. Nsalarm	two button warning station																
6. Triple pull	pull cord warning station																
7. External	use of external IIC device																
8. Syncalarm	synchronized two button warning station to be used with wireless Indicator																
Connections	<p>Built in pull cord switch on channel 0</p> <p>Built in reed switch in channel 2</p> <p>2 terminals for external emergency button</p> <p>2 terminals for extra external button</p> <p>2 terminals for input or output of 5V power supply</p> <p>2 terminals for external I2C bus</p> <p>2 sets of 2 terminals for external indicators</p> <p>4 terminals for auxiliary inputs/outputs</p> <p>4 UART pins, gnd, 5V, rxd, txd</p>																
Jumpers	<p>J1: choice internal external indicator 1</p> <p>J2: choice internal external indicator 2</p> <p>J13: prg, to set the NTM in wired programming mode. When removed programming is only possible via remote programming</p>																
Parameters	See application note 5; Manual for the nurse call binary station																
Indicators	LD1 Green LED																

NURSE CALL BINARY STATION

NCB	
	LD2 Red LED
Power supply and current consumption	Battery clip for CR123 battery Supply voltage 2.5 – 5 V Input/output power terminal: 5V Reverse polarity protected
NTM mode	The NCB mode is activated using function 7
Mounting instructions	For the best performance of the radio transceiver mount the PCB with the antenna upright when possible
Additional information	Datasheet NTM_3 Application note 1; Programming the NTM Application note 2; Ninthway high secure radio network Application note 5; Manual for the nurse call binary station